



Product data sheet REF 88.08.04

11.12.2025

| General data | | | |
|--|--|-------------------------------------|--|
| article number | 88.08.04 | bipol non-s bay fcps 20cm 1.3mm str | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 47848 | GTIN | 04046826423711 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 11-502 |
| Classification | 9 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,000 | Operating instructions | G880702 |
| Latex free | yes | Preparation instructions | MPK10 |
|  | | | |
| Medicon media | 450.00.37 (catalogue ENT) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.06.28 (brochure bipolar non-sticking fcps) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 20cm 1.3mm | Coating | powdered thermoplastic polyamide coating in various colors |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Others | SILBER_800 (AgCu_80-20) | | |