





Product data sheet REF 88.80.16

11.12.2025

| General data | | | |
|---|-----------------------------------|---------------------------|----------------|
| article number | 88.80.16 | shaft tube f bipolar fcps | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 47848 | GTIN | 04046826163914 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 11-502 |
| Classification | 9 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,037 | Operating instructions | G888010 |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | PPSU (Polyphenylsulfon) | | |