



Product data sheet REF 88.80.38

11.12.2025

General data			
article number	88.80.38	bipolar coag fcps rippled	
Unit	piece	Customs code (HS code)	90189084
GMDN	32684	GTIN	04046826307950
Risk category	2b (EU MDD)	UMDNS-Code	11-502
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,254	Operating instructions	G888038
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size		Coating	powdered thermoplastic polyamide coating in various colors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		
Handle	PEEK (Polyetheretherketon)		
Others	1.4301 (X5CrNi18-10)		