





Product data sheet REF 88.80.44

11.12.2025

General data			
article number	88.80.44	electr w rippled tip 34cm f bipolar coag fcps	
Unit	piece	Customs code (HS code)	90189084
GMDN	47848	GTIN	04046826308018
Risk category	2b (EU MDD)	UMDNS-Code	11-502
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,069	Operating instructions	G888038
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size	34cm	Coating	powdered thermoplastic polyamide coating in various colors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		
Others	PEEK (Polyetheretherketon)		