







Product data sheet REF 88.90.25

11.12.2025

General data			
article number	88.90.25	bipolator ang 18.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	47848	GTIN	04046826164096
Risk category	2b (EU MDD)	UMDNS-Code	15-579
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,015	Operating instructions	G870002
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
  			
size	18.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		