






## Product data sheet REF 91.05.00

11.12.2025

General data			
article number	91.05.00	endoscope Ø 5.0mm 300mm 0°	
Unit	piece	Customs code (HS code)	90189020
GMDN	12291	GTIN	04046826164935
Risk category	2a (EU MDD)	UMDNS-Code	12-291
Classification	5 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,325	Operating instructions	G776000
Latex free	yes	Preparation instructions	
			
Medicon media			
Complete picture		Detail images	
 			
size	Ø 5.0mm 300mm 0°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
<b>Main part</b>	1.4301 (X5CrNi18-10)		
<b>Screw</b>	1.4305 (X8CrNiS18-9)		
<b>Handle</b>	PEEK (Polyetheretherketon )		
<b>Others</b>	SAPHIRGLAS (Glas aus Saphirkristallen)		