







## Product data sheet REF 91.06.36

11.12.2025

General data			
article number	91.06.36	dilation set from 5mm to 11mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	37142	GTIN	04046826308261
Risk category	1r (EU MDR)	UMDNS-Code	11-254
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;01;aj;4W
Net weight (kg)	0,114	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
 			
size	5mm 11mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	SILIKON_MF670 (Silikon MF 670 70° Shore A)		