





Product data sheet REF 91.10.30

11.12.2025

General data			
article number	91.10.30	endoscope Ø 10mm 310mm 30°	
Unit	piece	Customs code (HS code)	90189020
GMDN	12291	GTIN	04046826165239
Risk category	2a (EU MDD)	UMDNS-Code	12-291
Classification	5 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,290	Operating instructions	G776000
Latex free	yes	Preparation instructions	
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 10mm 310mm 30°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	PEEK (Polyetheretherketon)		
Others	SAPHIRGLAS (Glas aus Saphirkristallen)		