



Product data sheet REF 91.14.12

11.12.2025

General data			
article number	91.14.12	automatic trocar sleeve Ø 12.5mm 98mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	37144	GTIN	04046826165581
Risk category	2a (EU MDD)	UMDNS-Code	14-154
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,136	Operating instructions	G911310
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 12.5mm 98mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		