



## Product data sheet REF 91.14.55

11.12.2025

General data			
article number	91.14.55	automatic trocar sleeve Ø 5.5mm 95mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	37144	GTIN	04046826165673
Risk category	2a (EU MDD)	UMDNS-Code	14-154
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,101	Operating instructions	G911310
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 5.5mm 95mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		