













Product data sheet REF 91.70.03

11.12.2025

General data			
article number	91.70.03	set of tubes single-use sterile p/10 f LS-i20+LS-i30	
Unit	box/10	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826167202
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,904	Operating instructions	
Latex free		Preparation instructions	KEIN MP
  			
Medicon media	450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
Storage conditions			
 60 °C	 90 %	 700 hPa	 1060 hPa
			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		