



Product data sheet REF 91.70.05

11.12.2025

General data			
article number	91.70.05	Clip set for motor cables, P/10 clips, sterilizable	
Unit	box/10	Customs code (HS code)	90189084
GMDN	46102	GTIN	04046826537630
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	KEIN MP
 			
Medicon media			
Complete picture		Detail images	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		