











Product data sheet REF 91.70.30

11.12.2025

| General data | | | |
|---|---|--|---|
| article number | 91.70.30 | infiltration pump LS-i30 cpl | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 60909 | GTIN | 04046826439989 |
| Risk category | N.A. (EU MDD) | UMDNS-Code | |
| Classification | NA (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 4,552 | Operating instructions | G917030 |
| Latex free | | Preparation instructions | |
|  | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
|  60 °C |  90 % |  1060 hPa |  |
| |  700 hPa |  | |
| Materials | | | |
| | Material number (DIN designation) | | |
| | | | |