



Product data sheet REF 93.00.04

11.12.2025

General data			
article number	93.00.04	wall tube VipAir w Synthes connector	
Unit	piece	Customs code (HS code)	90189084
GMDN	44505	GTIN	04046826168063
Risk category	2a (EU MDD)	UMDNS-Code	16-868
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	1,737	Operating instructions	G930000
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media			
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		