



Product data sheet REF 93.00.72

11.12.2025

General data			
article number	93.00.72	irrig tube 55mm f medium hdl	
Unit	piece	Customs code (HS code)	90189084
GMDN	44505	GTIN	04046826168254
Risk category	2a (EU MDD)	UMDNS-Code	16-868
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,001	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media			
Complete picture		Detail images	
			
size	55mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		