







## Product data sheet REF 93.00.85

11.12.2025

General data			
article number	93.00.85	double container 502x239x102mm	
Unit	piece	Customs code (HS code)	76169990
GMDN	12143	GTIN	04046826168315
Risk category	1 (EU MDR)	UMDNS-Code	13-730
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;04;XB
Net weight (kg)	5,052	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
 			
size	502x239x102mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	PPSU (Polyphenylsulfon)		