



Product data sheet REF 93.03.17

11.12.2025

General data			
article number	93.03.17	rosen burr L Ø 4.0mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	36249	GTIN	04046826168742
Risk category	2a (EU MDD)	UMDNS-Code	15-204
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,001	Operating instructions	G930110
Latex free	yes	Preparation instructions	MPK02
STERILE R			
Medicon media			
Complete picture		Detail images	
			
size	Ø 4.0mm	Coating	
Storage conditions			
5 °C  30 °C  10 % 60 %   			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		