



## Product data sheet REF 93.04.09

11.12.2025

General data					
article number	93.04.09	diamond dust burr XXL Ø 4.0mm			
Unit	piece	Customs code (HS code)	90189084		
GMDN	36249	GTIN	04046826168865		
Risk category	2a (EU MDD)	UMDNS-Code	15-204		
Classification	6 (EU MDD)	Base UDI-DI			
Net weight (kg)	0,001	Operating instructions	G930110		
Latex free	yes	Preparation instructions	MPK02		
<table border="1" style="margin: auto;"> <tr> <td style="padding: 5px;"><b>STERILE</b></td> <td style="padding: 5px;"><b>R</b></td> </tr> </table>				<b>STERILE</b>	<b>R</b>
<b>STERILE</b>	<b>R</b>				
Medicon media					
Complete picture		Detail images			
					
size	Ø 4.0mm	Coating	Artificially produced diamond coating on working ends of instruments		
Storage conditions					
<div style="display: flex; align-items: center; justify-content: space-around;"> <div style="text-align: center;">  <p>5 °C</p> </div> <div style="text-align: center;"> <p>30 °C</p> </div> <div style="text-align: center;">  <p>10 %</p> </div> <div style="text-align: center;"> <p>60 %</p> </div> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> </div>					
Materials					
	Material number (DIN designation)				
Main part	1.4112 (X90CrMoV18)				