

Ficha técnica REF 45.30.26

25.08.2026

Datos generales			
Referencia	45.30.26	repuesto tubo Satin	
Unidad	unid.	Número de aduana	90189084
GMDN	60909	GTIN	04046826534974
Clasificación	N.A. (EU MDR)	UMDNS-Code	
Clasificación	NA (EU MDR)	Base UDI-DI	
Peso neto (kg)	0,000	Manual de uso	
Latex free	si	Instrucciones de preparación	



Medicon medios

Imagen completa	Imagen detallada
-----------------	------------------

Date	Time	Location	Weather	Notes
10/10/2023	08:00	Forest	Clear	Initial observation
10/10/2023	08:15	Forest	Clear	Observation 1
10/10/2023	08:30	Forest	Clear	Observation 2
10/10/2023	08:45	Forest	Clear	Observation 3
10/10/2023	09:00	Forest	Clear	Observation 4
10/10/2023	09:15	Forest	Clear	Observation 5
10/10/2023	09:30	Forest	Clear	Observation 6
10/10/2023	09:45	Forest	Clear	Observation 7
10/10/2023	10:00	Forest	Clear	Observation 8
10/10/2023	10:15	Forest	Clear	Observation 9
10/10/2023	10:30	Forest	Clear	Observation 10
10/10/2023	10:45	Forest	Clear	Observation 11
10/10/2023	11:00	Forest	Clear	Observation 12
10/10/2023	11:15	Forest	Clear	Observation 13
10/10/2023	11:30	Forest	Clear	Observation 14
10/10/2023	11:45	Forest	Clear	Observation 15
10/10/2023	12:00	Forest	Clear	Observation 16
10/10/2023	12:15	Forest	Clear	Observation 17
10/10/2023	12:30	Forest	Clear	Observation 18
10/10/2023	12:45	Forest	Clear	Observation 19
10/10/2023	13:00	Forest	Clear	Observation 20
10/10/2023	13:15	Forest	Clear	Observation 21
10/10/2023	13:30	Forest	Clear	Observation 22
10/10/2023	13:45	Forest	Clear	Observation 23
10/10/2023	14:00	Forest	Clear	Observation 24
10/10/2023	14:15	Forest	Clear	Observation 25
10/10/2023	14:30	Forest		

Condiciones de almacenamiento

No se requieren condiciones especiales de conservación

Material

	N° de material (Denominación DIN)